

# ARTEch Laboratory

Amputee Restoration Technologies  
Midlothian, Texas 76065

309 W. Avenue F

Office 1-888-775-5501

Fax 1-972-775-2000

## Work Order

### Patient Profile

Prosthetic Company \_\_\_\_\_

Proth. Phone \_\_\_\_\_

Practitioner \_\_\_\_\_

Notes and Special Requirements \_\_\_\_\_

Patient Name \_\_\_\_\_

Description of amputation \_\_\_\_\_

Description of prosthesis(ae myo) \_\_\_\_\_

P.O. # \_\_\_\_\_

Practitioner' Signature \_\_\_\_\_ Date \_\_\_\_\_

***\*Mark sensitive areas on model & diagram***




Please complete all applicable profile/ diagram specifications and send this work order along with required models, photos, special instructions, etc.

#### Nose Colors

1. \_\_\_\_\_ Bridge of Nose
2. \_\_\_\_\_ Left Cheek
3. \_\_\_\_\_ Right Cheek
4. \_\_\_\_\_ Top Lip

#### Nose Measurements

5. \_\_\_\_\_ Cheekbone to Cheekbone
6. \_\_\_\_\_ Upper Lip to Corner of Right Eye
7. \_\_\_\_\_ Upper Lip to Corner of Left Eye
8. \_\_\_\_\_ Tip of Nose to Face
9. \_\_\_\_\_ Eye Corner to Eye Corner



**ARTech Laboratory, Inc.  
Checklist**

**Nose**

- Please read “Nose Casting Procedure”**
- Please read “Photo & Model Policy”**
- Photos of patient before amputation, if available. If not, photos of a family member or other person’s nose they would like.
- Photos taken ***before casting.***
- Dental stone model of residual
- Measurements completed on **Work Order**
- Color choices completed on **Work Order**

**ARTech Lab no longer fabricates nasal prostheses that require magnets or clips for retention. We only fabricate nasal prostheses with adhesive retention.**

**If patient has had a prosthesis previously, it would be helpful to us to have it, or at least a photo of it.** They usually want the new prosthesis to look similar to the one they are using now unless there is some problem with the current one.

**We do not need a lot of photos that are different colors. One or two prints of each angle, if they are color-true, will be sufficient.**

**It is not necessary to return the color chart with the models, unless you used one printed off our web site.**

**The very best way to ship stone models is to wrap them (after drying) in a bubble wrap.**