

ARTech Laboratory, Inc.
Checklist

Ear

- _____ **Please read “Photographic Policy”**
- _____ **Photos taken of both sides **before casting.****
- _____ **Please read “Ear Casting Procedure”**
- _____ **Lab Plaster model of sound side.**
- _____ **Dental stone model of residual affected side.**
- _____ **Color choices filled out on Work Order**

If patient has had a prosthesis previously, it would be helpful to us to have it, or at least a photo of it. They usually want the new prosthesis to look similar to the one they are using now unless there is some problem with the current one.

We do not need a lot of photos that are different colors. One or two prints of each angle, if they are color-true, will be sufficient.

It is not necessary to return the color chart with the models, unless it was printed off our web site.