

ARTech Laboratory, Inc.
Checklist

Ear

- _____ Please read “**Ear Casting Procedure**”
- _____ **Please read “Photographic Policy”**
- _____ Photos taken of both sound & affected sides **before casting.**
- _____ Lab plaster model of sound side.
- _____ Dental stone model of residual affected side.
- _____ Color choices filled out on Work Order
- _____ Measurements filled out on Work Order

ARTech Lab no longer fabricates ear prostheses that require magnets or clips for retention. We only fabricate ear prostheses that have adhesive retention.

If patient has had a prosthesis previously, it would be helpful to us to have it, or at least a photo of it. They usually want the new prosthesis to look similar to the one they are using now unless there is some problem with the current one.

We do not need a lot of photos that are different colors. One or two prints of each angle, if they are color-true, will be sufficient.

It is not necessary to return the color chart with the models, unless it was printed off our web site.